

Disaster/Transport Bag and ID Kit

Busse has added a Disaster Bag and a Disaster Bag Identification Kit to our existing line of post mortem products. This Disaster Bag is made to meet Government specifications as well as OSHA safety standards.

Catalog No. 920

Packed 5 ea./carton

34" x 96" Black Disaster Bag with ID Kit

- Meets US Government Specifications and Department of Defense (DOD) Specifications NSN 9930-01-331-6244
- Complies with OSHA Regulation 3130, Occupational Exposure To Bloodborne Pathogens — Precautions for Emergency Responders
- Heat sealed seams prevent leakage of fluid from the sides during transport
- Envelope style easy access zipper opening
- 20 mil black scrim reinforced PVC for strength
- 6 handles for ease of carrying
- Personal effects bag
- Latex free



Catalog No. 921

Packed 20 ea./carton

Body Bag Identification Kit

- Kit comes in clear plastic pouch with self adhesive strips
- Includes 1 red biohazard tag, 3 yellow contamination tags, 3 identification tags and self adhesive document pouch with biohazard symbol
- Maintain identification and chain of custody
- Latex free



busse

Hospital Disposables

P.O. Box 11067, Hauppauge, NY 11788-0920

Tel: 800.654.6526 • Fax: 631.435.4721

www.busseinc.com • mail@busseinc.com

© 2010 Busse Hospital Disposables

E.O.E.—M/F/D/V

920 • Product specifications subject to change without notice.