

STANDARD PAIN MANAGEMENT TRAYS

TRAY COMPONENTS		Basic Pain Trays		Epidural Single Shot				Epidural Continuous	Nerve Block			
		674	675	670	677	678	679	682	671	672	733	734
KEY COMPONENTS	18g x 3.5" Tuohy Epidural Needle, Fixed Wings, Metal Stylet						●	●				
	18g x 3.5" Tuohy Epidural Needle, Fixed Wings			●				●				
	20g x 3.5" Tuohy Epidural Needle, Fixed Wings, Metal Stylet				●	●						
	22g x 3.5" Quincke Spinal Needle									●		
	25g x 3.5" Quincke Spinal Needle											●
	10mL Busse LOR Syringe, Luer Lock Tip	●		●		●						
	10mL Busse LOR Syringe, Luer Slip Tip				●		●					
	5mL Glass LOR Syringe, Metal Luer Slip Tip			●				●				
	20g Closed-End Epidural Catheter, Thread Assist Guide, Catheter Pad, Label and Tuohy Borst Adapter							●				
	Epidural Filter							●				
PHARMACEUTICALS	5mL, 1% Lidocaine Ampule		●	●				●		●		
	5mL, 1.5% Lidocaine with Epinephrine Ampule							●				
	1mL, 0.1% Epinephrine Ampule							●				
	10mL, 0.9% Sodium Chloride Ampule		●	●				●				
PREP COMPONENTS	PVP Solution	●	●	●	●	●	●	●	●	●	●	●
	Sponge Sticks	3	3	3	3	3	3	3	3	3	3	3
	Absorbent Towel	●	●	●	●	●	●	●	●	2	●	●
	Blue Fenestrated Drape with Adhesive Strip	●	●		●	●	●	●				
	Clear Fenestrated Drape with Adhesive Strip			●				●		●	●	●
	Skin Marker									●		
	Flexible Ruler									●		
	Gauze Dressings	4	4	4	4	4	4	4	4	4	4	4
PROCEDURAL COMPONENTS	1mL Luer Lock Syringe									●		
	3mL Luer Lock Syringe	2	2	●	2	2	2	2	●	●	2	2
	5mL Luer Lock Syringe	●	●		●	●	●	●		●	●	●
	10mL Luer Lock Syringe	●	●		●	●	●	●			●	●
	20mL Luer Lock Syringe			●				●				
	20mL Control Syringe									●		
	18g x 1.5" Needle	●	●	●	●	●	●	●	●	●	●	●
	25g x 1.5" Needle	●	●	●	●	●	●	●	●	●	●	●
	27g x 1.25" Needle	●	●		●	●	●	●			●	●
	18g x 1.5" Filter Needle		●	●				●		●		
	Needle Stick Block			●				●		●		
	Busse Saf-T-Cube™	●	●		●	●	●				●	●
	Adhesive Bandage			●				●		●		
Medication Labels	●	●		●	●	●	●			●	●	

UNIVERSAL PREP TRAY		CATALOG NO.	730
TRAY COMPONENTS	PVP Solution	●	
	Sponge Sticks	3	
	Absorbent Towel	●	
	Blue Fenestrated Drape with Adhesive Strip	●	
	Gauze Dressings	4	
	Single Compartment Tray	●	

Components listed above are packaged in pre-formed prep and multi-compartment trays and then wrapped in CSR overwrap. All packed 10 per case.

busse
Hospital Disposables

75 Arkay Drive • Hauppauge, NY 11788-0920 • Tel: 800-645-6526
Fax: 631-435-4721 • www.busseinc.com • mail@busseinc.com

© 2013 Busse Hospital Disposables
Product specifications are subject to change.
E.O.E. - M/F/D/V
Standard Pain Management Trays



The universal prep tray is wrapped in a CSR overwrap.

Packed 10 per case.

STANDARD PAIN MANAGEMENT TRAYS



**BUSSE KNOWS
WHAT'S ON YOUR
WISH LIST**

...Standard trays with the most requested components.



busse
Hospital Disposables

Busse is proud to be a certified woman owned and operated company.

