

PREMIUM QUALITY Single-Use Instruments

Featuring our new
SUTURE-GLIDE™
Webster
Needle Holder

NO SNAGS!
NO HOLD UPS!
NO HASSLES!



- ER Laceration Tray with SUTURE-GLIDE™ Webster Needle Holder. Catalog No. 768. Packed sterile 20/case.
- SUTURE-GLIDE™ Webster Needle Holder also available individually packaged. Catalog No. 781. Packed sterile 25/case.

busse
Hospital Disposables


CUSTOMIZE YOUR ER LACERATION TRAYS OR SUTURE REMOVAL KITS.

Include any of our new premium quality, single-use instruments.

The traditional "box-lock" hinge design.



This is just a sampling of our new premium quality, single-use instruments.



Unlike the basic "box-lock" hinge, the uniquely designed SUTURE-GLIDE™ hinge features smooth beveled edges that allow the looped suture to gently and smoothly slide over the instrument jaws without snagging or tearing.

Exceptional design.

Remarkable value.

- Uniform smooth edges and bevels highlight the symmetrical design.
- Larger finger loops.
- Satin finish.

Busse offers a full line of premium quality, single-use instruments such as:

- Needle Holders
- Scissors
- Forceps
- Hemostats

Please go to our website or scan the QR code below on your mobile device to see our full line of ER Laceration Trays and Suture Removal Kits.




Sharp/Blunt Scissors



Iris Scissors



SUTURE-GLIDE™
Webster
Needle Holder
Smooth



Adson Tissue Forceps



Adson
Serrated Forceps



busse

Hospital Disposables

P.O. Box 11067 • Hauppauge, NY 11788-0920 • Tel: 800-645-6526 • Fax: 631-435-4721
www.busseinc.com • mail@busseinc.com

For more details, contact our custom kit department at 800-645-6526 or ask your Busse sales representative.

768 / Busse Hospital Disposables / Product specifications are subject to change.