



OEM SURGICAL DRAPES FOR KIT PACKERS

HEALTHCARE  
HAS A NEED TO  
CUSTOMIZE.

BUSSE HAS  
A CUSTOMIZED  
SOLUTION.

**busse**  
Hospital Disposables

# CUSTOM QUALITY IS STANDARD AT BUSSE.

At Busse, we are kit packers ourselves, so we fully understand the production issues you face on a daily basis. Busse has been manufacturing drapes for over 30 years, and our focus is clear. We process orders quickly and accurately... expedite shipment turnaround... and maintain exceptional quality assurance levels.

We are pleased to offer a wide variety of both standard and custom-configured procedural drapes.

- **Hard-to-find drapes are now available.**
- **A variety of drape materials are offered in absorbent, non-absorbent, clear and translucent options.**
- **Choose from a broad range of fenestration shapes and sizes. Fenestrations can be customized with and without adhesive borders or strips.**
- **Stringent quality control assures that every drape meets your requirements.**
- **Drapes are manufactured domestically in our manufacturing facility— your assurance of quality and prompt delivery.**



THE LISTING BELOW IS REPRESENTATIVE OF THE MANY STYLES OF DRAPES THAT WE CAN SUPPLY.

IF YOU DON'T SEE WHAT YOU NEED, PLEASE CONTACT OUR OEM DEPARTMENT AT 800-645-6526 FOR A CUSTOMIZED SOLUTION.

MINOR PROCEDURE / GENERAL SURGERY / OR

	CAT#	MATERIAL	SIZE (INCHES)	COLOR	DESCRIP.	PERFORMANCE CHARACTERISTIC/ BARRIER LEVEL	PKG. (CASE)
	6830	Plastic	17 x 11	Matte Finish	Drape with adhesive strip/fan-folded	Impervious L4	200
	6840	Plastic	17 x 23	Matte Finish	Drape with adhesive strip/fan-folded	Impervious L4	200
	6850	Plastic	15 x 15	Matte Finish	Drape with 2 1/2" fenestration & adhesive patch	Impervious L4	200
	6855	Plastic	15 x 15	Matte Finish	Drape with 2 1/2" fenestration and slit	Impervious L4	200
	6880	Clear Plastic	47 x 51	Clear	"U" Drape with Matte Finish	Impervious L4	200
	6711	Reinforced Tissue	19 1/2 x 23	Blue	Nylon Scrim Reinforced Towel	Absorbent N/A	500
	6935	Reinforced Tissue	19 1/2 x 28	Blue	Nylon Scrim Reinforced Towel	Absorbent N/A	500
	6861	Laminated Nonwoven	22 x 26	Blue	Laminated Nonwoven Drape	Absorbent/Impervious L4	500
	6870	Laminated Nonwoven	22 x 26	Blue	Laminated Nonwoven Drape with 2 1/2" fenestration & adhesive patch	Absorbent/Impervious L4	500
	6808	SMS	15 1/2 x 25	Blue	Utility Drape	Repellent L3	500

\*T = Tissue  
 \*PE = Polyethylene  
 \*PP = Polypropylene

Can be sterilized by EO














MINOR PROCEDURE / GENERAL SURGERY / OR (CONTINUED)

	CAT#	MATERIAL	SIZE (INCHES)	COLOR	DESCRIP.	PERFORMANCE CHARACTERISTIC/ BARRIER LEVEL	PKG. (CASE)
	6820	SMS	15 1/2 x 25	Blue	Utility Drape with adhesive strip	Repellent L3	500
	6825	SMS	22 x 26	Blue	Drape with 4" x 4" square fenestration & adhesive patch	Repellent L3	250
	695	*T/PP/T	18 x 26	White	Drape with Blue PP Film	Absorb/Imperv/Absorb L4	500
	694	*T/PP/T	18 x 26	White	Drape with Blue PP Film & 3" fenestration	Absorb/Imperv/Absorb L4	300
	6960	*T/PE/T	18 x 26	White	Towel/Drape	Absorb/Imperv/Absorb L4	500
	6970	*T/PE/T	18 x 26	White	Drape with 3" fenestration	Absorb/Imperv/Absorb L4	500
	6598	*T/PE/T	18 x 26	White	Drape with 3" fenestration & adhesive patch	Absorb/Imperv/Absorb L4	500
	6600	*T/PE/T	24 x 36	White	Drape with 4" fenestration & adhesive patch	Absorb/Imperv/Absorb L4	250

OPHTHALMOLOGY

	6850	Plastic	15 x 15	Matte Finish	Drape with 2 1/2" fenestration & adhesive patch	Impervious L4	200
	6890	Plastic	51 x 51	Matte Finish	Drape with 3 7/8" x 4 7/8" rectangular incise area & fluid collection pouch	Impervious L4	200

## PAIN MANAGEMENT

	CAT#	MATERIAL	SIZE (INCHES)	COLOR	DESCRIP.	PERFORMANCE CHARACTERISTIC/ BARRIER LEVEL	PKG. (CASE)
	6715	Plastic	24 x 30	Matte Finish	Drape with 4" fenestration & tape strip	Impervious L4	250
	6741	Clear Plastic	21 x 30	Clear	Drape with 4" fenestration & adhesive patch	Impervious L4	250
	6740	Clear Plastic	21 x 36	Clear	Drape with 4" fenestration & adhesive patch	Impervious L4	250
	6720	Clear Plastic	24 x 30	Clear	Drape with 4" fenestration & tape strip	Impervious L4	250
	6735	Clear Plastic	24 x 30	Clear	Drape with 5" fenestration & tape strip, banded	Impervious L4	250
	6745	Clear Plastic	24 x 30	Clear	Drape with 4"x 8" rectangular fenestration & tape strip	Impervious L4	250
	6750	Clear Plastic	24 x 30	Clear	Drape with 4" fenestration, adhesive patch & 2 adhesive tabs	Impervious L4	250
	6760	Clear Plastic	24 x 36	Clear	Drape with 4" fenestration, adhesive patch & 2 adhesive tabs	Impervious L4	250
	6730	Clear Plastic	28 x 30	Clear	Drape with 5" fenestration & tape strip	Impervious L4	250
	6725	Nonwoven	24 x 30	Blue	Drape with 4" fenestration & tape strip	Absorbent/ Repellent L3	250
	6825	SMS	22 x 26	Blue	Drape with 4"x 4" square fenestration & adhesive patch	Repellent L3	250

### COLOR KEY

### BARRIER PERFORMANCE LEVELS

IMPERVIOUS	4
ABSORBENT/IMPERVIOUS	4
ABSORB./IMPERV./ABSORB.	4
REPELLENT	3
ABSORBENT/REPELLENT	3
ABSORBENT	N/A
ADHESIVE AREA	N/A

Busse is ISO Certified

Level of barrier performance is based on the AAMI PB 70\*

All Busse drapes meet flammability standard 16 CFR Park 1610 – Type A, Class 1\*

Custom folding and banding is available

All Busse drapes are latex free

All Busse drapes are EO sterilizable

Customized drapes may require a minimum quantity. Please contact our OEM department at 800-645-6526.

\*Test results are available upon request for performance and flammability

**busse**  
Hospital Disposables

## CUSTOM DESIGN. CUSTOM FABRICATION. CONSISTENT QUALITY.

Drape size, configuration, fenestration variety and adhesive application have always been limited for kit packers until now. Busse offers greater options to meet your customers' requirements.

- **Fenestrations are die-cut for a smooth, clean edge.**
- **Adhesive applications are precisely monitored to assure trouble-free unfolding and perfect drape placement.**
- **Our manufacturing process can adjust fenestration location and adhesive placement based on specific requirements.**



Traditional circular fenestrations measuring from 5/8 inches to over 5 inches in diameter are available.




Die-cut square fenestrations assure clean border edges, and can be configured up to 6 inches.



Large rectangular fenestrations can now be custom cut to suit a greater operating field, thus eliminating the need to reposition the drape.

Busse drapes are manufactured in a variety of quality materials and layer configurations. All of our surgical drapes are fabricated with the utmost care, assuring you and your customers of the ultimate in protection and ease of use.



Busse's advanced manufacturing process provides for the precise placement of adhesive patches, strips and tabs.





Hospital Disposables

P.O. Box 11067 • Hauppauge, NY 11788-0920 • Tel: 800-645-6526 • Fax: 631-435-4721  
[www.busseinc.com](http://www.busseinc.com) • [mail@busseinc.com](mailto:mail@busseinc.com)

© 2009 Busse Hospital Disposables  
Product specifications are subject to change.  
E.O.E. - M/F/D/V  
OEM Drapes